ProphyPal®

Hygiene Handpiece



Sure-fit™ Nosecone

Extended nosecone provides excellent stability for all disposable prophy angles.

360° Swivel and Grip Nosecone

Allows for maximum control and easy transfer between buccal and lingual surfaces.

Removable Nosecone

Easily removable nosecone allows for efficient sterilization. Fully autoclave your nosecone and handpiece up to 275°F/135°C.

Power and Torque

0 - 5,000 RPM motor operates at optimal speed and torque at 40-50 PSI. The whisper quiet air powered operation reduces patient and operator stress.

Ergonomic

The lightweight and balanced construction provides excellent ergonomic positions for all hand sizes.

Available in Classic Silver, Prophy Pink or Wild Thing Green Colors



